

Small Fragment

Small Compression Plates 3.5 mm

Material: Stainless Steel
 Screw diameter: 3.5 mm
 Plate thickness: 3.0 mm
 Plate width: 10.0 mm
 Hole spacing: 12.0 mm / 16.0 mm
 Application: usually used for radius, ulna, distal humerus and fibula fixation, and forearm osteotomies.

Fabricated in Germany by Ortho Select® GmbH Implant Technology, ISO 13485:2003 Certified.



Not to scale

Holes	Length (mm)	Part No.	Length (mm)	Minimal Contact Part No.
2	25	6010-0-4002		
3	37	6010-0-4003		
4	49	6010-0-4004	51	6011-0-4004
5	61	6010-0-4005	64	6011-0-4005
6	73	6010-0-4006	77	6011-0-4006
7	85	6010-0-4007	90	6011-0-4007
8	97	6010-0-4008	103	6011-0-4008
9	109	6010-0-4009	116	6011-0-4009
10	121	6010-0-4010	129	6011-0-4010
11	133	6010-0-4011	142	6011-0-4011
12	145	6010-0-4012	155	6011-0-4012
14	169	6010-0-4014		
16	194	6010-0-4016		
18	217	6010-0-4018		