



Small Fragment

One-Third Tubular Plates

Material: Stainless Steel
 Screw diameter: 3.5 mm
 Plate thickness: 1.0 mm
 Plate width: 9.0 mm
 Hole spacing: 12.00 mm - 16.0 mm
 Application: usually used for distal ulna and fibula fixation.

Fabricated in Germany by Ortho Select® GmbH Implant Technology, ISO 13485:2003 Certified.



Not to scale

Holes	Length mm	Without Collar Part No.	With Collar Part No.
2	25	6022-0-4002	6022-2-4002
3	37	6022-0-4003	6022-2-4003
4	49	6022-0-4004	6022-2-4004
5	61	6022-0-4005	6022-2-4005
6	73	6022-0-4006	6022-2-4006
7	85	6022-0-4007	6022-2-4007
8	97	6022-0-4008	6022-2-4008
9	109	6022-0-4009	6022-2-4009
10	121	6022-0-4010	6022-2-4010
11	133	6022-0-4011	6022-2-4011
12	145	6022-0-4012	6022-2-4012