

Small Fragment

Small T-Plates

Material: Stainless Steel
 Screw diameter: 3.5 mm
 Plate thickness: 1.2 mm
 Plate width: 11.0 mm
 Application: usually used for distal radius fixation and osteotomies.

Fabricated in Germany by Ortho Select[®] GmbH Implant Technology, ISO 13485:2003 Certified.



Shaft Holes	Head Holes	Length (mm)	Oblique Part No.
3	3	52	6052-0-4033
4	3	63	6052-0-4043
5	3	74	6052-0-4053



Shaft Holes	Head Holes	Length (mm)	90° Part No.
3	3	50	6050-0-4033
4	4	56	6050-0-4044
5	3	67	6050-0-4053
6	4	78	6050-0-4064

Application: usually used for distal radius and tibial plateau fixation, and distal radius osteotomies.