



Small Fragment

Cloverleaf Plates

Material: Stainless Steel
Screw diameter: 3.5 mm
Plate thickness: 2.0 mm
Plate width: 15.0 mm
Application: usually used for distal tibial and proximal humerus fixation

Fabricated in Germany by Ortho Select® GmbH Implant Technology, ISO 13485:2003 Certified.



Not to scale

Holes	Length (mm)	Part No.
3	88	6054-0-4003
4	104	6054-0-4004
5	120	6054-0-4005
6	136	6054-0-4006