

Small Fragment

Cortex Screws 3.5 mm

Material:	Stainless Steel
Thread diameter:	3.5 mm
Core diameter:	2.4 mm
Head diameter:	6.0 mm
Hex diameter:	2.5 mm
Thread pitch:	1.25 mm
Pilot drill:	2.5 mm
Gliding hole drill:	3.5 mm
Tap:	3.5 mm
Screwdriver hex:	2.5 mm
Application:	usually used for interfragmentary compression according to the lag-technique, for plate fixation, or for both.

Fabricated in Germany by Ortho Select[®] GmbH Implant Technology; ISO 13485:2003 Certified.



Length (mm)	3.5 mm Standard Part No.	3.5 mm Self-Tapping Part No.
10	6080-0-4010	6081-0-4010
12	6080-0-4012	6081-0-4012
14	6080-0-4014	6081-0-4014
16	6080-0-4016	6081-0-4016
18	6080-0-4018	6081-0-4018
20	6080-0-4020	6081-0-4020
22	6080-0-4022	6081-0-4022
24	6080-0-4024	6081-0-4024
26	6080-0-4026	6081-0-4026
28	6080-0-4028	6081-0-4028
30	6080-0-4030	6081-0-4030
32	6080-0-4032	6081-0-4032
34	6080-0-4034	6081-0-4034
36	6080-0-4036	6081-0-4036

Small Fragment

Cortex Screws 3.5 mm



Length (mm)	3.5 mm Standard Part No.	3.5 mm Self-Tapping Part No.
38	6080-0-4038	6081-0-4038
40	6080-0-4040	6081-0-4040
45	6080-0-4045	6081-0-4045
50	6080-0-4050	6081-0-4050
55	6080-0-4055	6081-0-4055
60	6080-0-4060	6081-0-4060
65	6080-0-4065	6081-0-4065
70	6080-0-4070	6081-0-4070
75	6080-0-4075	6081-0-4075
80	6080-0-4080	6081-0-4080
85	6080-0-4085	6081-0-4085
90	6080-0-4090	6081-0-4090
95	6080-0-4095	6081-0-4095
100	6080-0-4100	6081-0-4100
105	6080-0-4105	6081-0-4105
110	6080-0-4110	6081-0-4110
115	6080-0-4115	6081-0-4115
120	6080-0-4120	6081-0-4120