

Small Fragment

Cancellous Screws


Material: Stainless Steel per ISO 5832-1

Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.



Thread diameter:	3.5 mm	4.0 mm	4.0 mm
Head diameter:	6.0 mm	6.0 mm	6.0 mm
Hex diameter:	2.5 mm	2.5 mm	2.5 mm
Thread pitch:	1.75 mm	1.75 mm	1.75 mm
Pilot drill:	2.5 mm	2.5 mm	2.5 mm
Tap:	3.5 mm	4.0 mm	4.0 mm
Screwdriver hex:	2.5 mm	2.5 mm	2.5 mm
	3.5 mm Standard Part No.	4.0 mm Fully Threaded Part No.	4.0 mm Partially Threaded Part No.
Length (mm)			
10	6083-0-4010	6084-0-5010	6085-0-5010
12	6083-0-4012	6084-0-5012	6085-0-5012
14	6083-0-4014	6084-0-5014	6085-0-5014
16	6083-0-4016	6084-0-5016	6085-0-5016
18	6083-0-4018	6084-0-5018	6085-0-5018
20	6083-0-4020	6084-0-5020	6085-0-5020
22	6083-0-4022	6084-0-5022	6085-0-5022
24	6083-0-4024	6084-0-5024	6085-0-5024
26	6083-0-4026	6084-0-5026	6085-0-5026
28	6083-0-4028	6084-0-5028	6085-0-5028
30	6083-0-4030	6084-0-5030	6085-0-5030

Continued on page 2



Small Fragment

Cancellous Screws



Length (mm)	3.5 mm Standard Part No.	4.0 mm Fully Threaded Part No.	4.0 mm Partially Threaded Part No.
32	6083-0-4032		6085-0-5032
34	6083-0-4034	6084-0-5034	
35		6084-0-5035	6085-0-5035
36	6083-0-4036		
38	6083-0-4038		
40	6083-0-4040	6084-0-5040	6085-0-5040
45	6083-0-4045	6084-0-5045	6085-0-5045
50	6083-0-4050	6084-0-5050	6085-0-5050
55	6083-0-4055	6084-0-5055	6085-0-5055
60	6083-0-4060	6084-0-5060	6085-0-5060