



## Mini Fragment

### T-Plate 1.5 mm

Material: Stainless Steel  
Screw diameter: 1.5 mm  
Shaft thickness: 0.9 mm  
Shaft width: 3.8 mm  
Hole spacing: 5.0 mm  
Application: usually used for metatarsal, metacarpal, and phalanx fixation.

*Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.*



Shaft Holes	Head Holes	Part No.
9	3	6009-0-1509
9	4	6009-0-1510