



Mini Fragment

H-Plate 1.5 mm

Material: Stainless Steel
Screw diameter: 1.5 mm
Shaft thickness: 0.9 mm
Outer dimensions : 9.0 mm x 8.2 mm
Hole spacing: 5.0 mm x 4.2 mm
Application: usually used for metatarsal and metacarpal fixation.

Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.



Not to scale

Holes	Part No.
4	6009-0-1530