



Mini Fragment

Condylar Plate 1.5 mm

Material:	Stainless Steel
Screw diameter:	1.5 mm
Shaft thickness:	0.9 mm
Shaft width:	3.5 mm
Hole spacing:	5.0 mm
Length:	26.0 mm
Application:	usually used for metacarpal osteotomies.

Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.



Holes		Part No.
6	Right	6009-0-1550
6	Left	6009-0-1551