



## Mini Fragment

### Bone Plate 2.0 mm

Material:	Stainless Steel
Screw diameter:	2.0 mm
Shaft thickness:	1.0 mm (4-6 holes)
Shaft thickness:	1.5 mm (7-8 holes)
Shaft width:	5.0 mm
Hole spacing:	5.0 mm
Application:	usually used for fixation of hand and foot.

*Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.*



Holes	Length	Part No.
4	22	6009-0-2004
5	27	6009-0-2005
6	32	6009-0-2006
7	37	6009-0-2007
8	42	2009-0-2008