



## Mini Fragment

### T-Plate 2.0 mm

Material: Stainless Steel  
Screw diameter: 2.0 mm  
Shaft thickness: 1.2 mm  
Shaft width: 5.0 mm  
Hole spacing: 5.0 mm  
Application: usually used for proximal phalanx fixation and osteotomies.

*Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.*



Shaft Holes	Head Holes	Part No.
9	3	6009-0-2016
9	4	6009-0-2017