



Mini Fragment

Straight Mini Plate 2.0 mm

Material: Stainless Steel
Screw diameter: 1.5 mm
Shaft thickness: 1.2 mm
Shaft width: 5.0 mm
Hole spacing: 6.0 mm
Application: usually used for phalanx fixation.

Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.



Holes	Length (mm)	Part No.
3	17	6009-0-2103
4	23	6009-0-2104
5	29	6009-0-2105
6	35	6009-0-2106