



Mini Fragment

L-Plate 2.0 mm

Material: Stainless Steel
Screw diameter: 1.5 mm
Shaft thickness: 1.0 mm
Shaft width: 5.0 mm
Hole spacing: 6.0 mm
Application: usually used for metacarpal, phalanx fixation, and osteotomies.

Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.



Holes		Part No.
2	Oblique - Right	6009-0-2002
2	90° - Right	6009-0-2012
2	Oblique - Left	6009-0-2003
2	90° - Left	6009-0-2011