



## Mini Fragment

### T-Plate 2.0 mm

Material:	Stainless Steel
Screw diameter:	2.0 mm
Shaft thickness:	1.0 mm
Shaft width:	5.0 mm
Hole spacing:	6.0 mm
Application:	usually used for metatarsal, metacarpal, phalanx fixation, and osteotomies.

*Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.*



*Not to scale*

Holes	Head Holes	Part No.
2	2	6009-0-2010