



Mini Fragment

H-Plate 2.0 mm

Material: Stainless Steel
Screw diameter: 2.0 mm
Shaft thickness: 0.9 mm
Outer dimensions : 11.0 mm x 12.0 mm
Hole spacing: 7.0 mm x 6.0 mm
Application: usually used for metatarsal and phalanx fixation.

Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.



Not to scale

Holes	Part No.
4	6009-0-2018