



Mini Fragment

Condylar Plate 2.0 mm

Material:	Stainless Steel
Screw diameter:	2.0 mm
Shaft thickness:	1.2 mm
Shaft width:	5.0 mm
Hole spacing:	5.0 mm
Length:	26.0 mm
Application:	usually used for metacarpal osteotomies.

Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.



Holes		Part No.
6	Right	6009-0-2050
6	Left	6009-0-2051