



Mini Fragment

Compression Plate 2.7 mm

Material:	Stainless Steel
Screw diameter:	2.7 mm
Shaft thickness:	2.0 mm (2-6 holes)
Shaft thickness:	2.5 mm (6-12 holes)
Shaft width:	8.0 mm
Hole spacing:	8.0 mm
Application:	usually used for metatarsal, metacarpal fixation, and osteotomies.

Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.



Holes	Length (mm)	Part No.
2	20	6009-0-2702
3	26	6009-0-2703
4	36	6009-0-2704
5	44	6009-0-2705
6	52	6009-0-2706
7	60	6009-0-2707
8	68	6009-0-2708
9	76	6009-0-2709
10	84	6009-0-2710
12	100	6009-0-2712