



Mini Fragment

L-Plate 2.7 mm

Material: Stainless Steel
Shaft thickness: 1.2 mm
Shaft width: 7.0 mm
Hole spacing: 8.0 mm
Application: usually used for metacarpal, phalanx fixation, and osteotomies.

Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.



Holes		Part No.
3	Oblique - Right	6009-0-2714
3	90° - Right	6009-0-2716
3	Oblique - Left	6009-0-2715
3	90° - Left	6009-0-2717