



Mini Fragment

Quarter Tubular Plates 2.7 mm

Material: Stainless Steel
Screw diameter: 2.7 mm
Shaft thickness: 1.0 mm
Shaft width: 7.0 mm
Hole spacing: 8.0 mm

Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.



Not to scale

Holes	Length (mm)	Part No.
3	23	6009-0-2733
4	31	6009-0-2734
5	39	6009-0-2735
6	47	6009-0-2736
7	55	6009-0-2737
8	63	6009-0-2738
10	71	6009-0-2740