



Mini Fragment

T-Plate 2.7 mm

Material: Stainless Steel
Screw diameter: 2.7 mm
Shaft thickness: 1.0 mm
Shaft width: 5.0 mm
Application: usually used for metatarsal, metacarpal, phalanx fixation, and osteotomies.

Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.



Not to scale

Holes	Length (mm)	Part No.
3	32	6009-0-2723
4	40	6009-0-2724
5	48	6009-0-2725
6	56	6009-0-2726
7	54	6009-0-2727
8	72	6009-0-2728