

Mini Fragment

Cortex Screws

Material: Stainless Steel

Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.



Not to scale, not in relative proportion.

Thread diameter:	1.5 mm	2.0 mm	2.0 mm
Core diameter:	1.0 mm	1.3 mm	1.3 mm
Thread pitch:	0.5 mm	0.6 mm	0.6 mm
Head diameter:	3.0 mm	4.0 mm	4.0 mm
Hex diameter:	1.5 mm	1.5 mm	1.5 mm
Pilot drill:	1.1 mm	1.5 mm	1.5 mm
Gliding hole drill:	1.5 mm	2.0 mm	2.0 mm
Tap:	1.5 mm	2.0 mm	2.0 mm
Screwdriver hex:	1.5 mm	1.5 mm	1.5 mm
Application:	Usually used to fix fracture fragments in or at phalanges, metacarpal, and small bones of the foot.	Usually used to fix fracture fragments in the proximal and middle phalanges, metacarpal, and small bones of the foot.	Usually used in the metacarpal and metatarsal bones. In some indications also used in radius and ulna.
Length (mm)	1.5 mm Part No.	2.0 mm Standard Part No.	2.0 mm Self-Tapping Part No.
6	6077-0-4006	6078-0-4006	
7	6077-0-4007		
8	6077-0-4008	6078-0-4008	
9	6077-0-4009		
10	6077-0-4010	6078-0-4010	6078-0-5010
11	6077-0-4011		

Mini Fragment

Cortex Screws



Not to scale, not in relative proportion.

Length (mm)	1.5 mm	2.0 mm Standard Part No.	2.0 mm Self-Tapping Part No.
12	6077-0-4012	6078-0-4012	6078-0-5012
14	6077-0-4014	6078-0-4014	6078-0-5014
16	6077-0-4016	6078-0-4016	6078-0-5016
18	6077-0-4018	6078-0-4018	6078-0-5018
20	6077-0-4020	6078-0-4020	6078-0-5020
22	6077-0-4022	6078-0-4022	6078-0-5022
24	6077-0-4024	6078-0-4024	6078-0-5024
26		6078-0-4026	6078-0-5026
28		6078-0-4028	6078-0-5028
30		6078-0-4030	6078-0-5030
32		6078-0-4032	6078-0-5032
34		6078-0-4034	6078-0-5034
36		6078-0-4036	6078-0-5036
38		6078-0-4038	6078-0-5038