



Mini Fragment

Cortex Screws 2.7 mm

Material:	Stainless Steel per ISO 5832-1
Thread diameter:	2.7 mm
Core diameter:	1.9 mm
Head diameter:	5.0 mm
Hex diameter:	2.5 mm
Pilot drill:	2.0 mm
Gliding hole drill:	2.7 mm
Tap:	2.7 mm
Screwdriver hex:	2.5 mm

Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.



Length (mm)	2.7 mm Standard Part No.	2.7 mm Self-Tapping Part No.
6	6079-0-4006	
8	6079-0-4008	
10	6079-0-4010	6079-0-5010
12	6079-0-4012	6079-0-5012
14	6079-0-4014	6079-0-5014
16	6079-0-4016	6079-0-5016
18	6079-0-4018	6079-0-5018
20	6079-0-4020	6079-0-5020
22	6079-0-4022	6079-0-5022
24	6079-0-4024	6079-0-5024
26	6079-0-4026	6079-0-5026
28	6079-0-4028	6079-0-5028
30	6079-0-4030	6079-0-5030
32	6079-0-4032	6079-0-5032
34	6079-0-4034	6079-0-5034
36	6079-0-4036	6079-0-5036
38	6079-0-4038	6079-0-5038
40	6079-0-4040	6079-0-5040