



## LARGE FRAGMENT

### T-Buttress Plates

Material:	Stainless Steel
Plate thickness:	2.0 mm
Plate width:	17.0 mm
Head width:	33.0 mm
Application:	usually used for lateral tibial plateau fixation

*Fabricated in Germany by Ortho Select® Implant Technology an ISO 13485:2003 Certified Company.*



Holes	Length (mm)	Part No.
4	80	6016-0-6004
5	96	6016-0-6005
6	112	6016-0-6006
8	144	6016-0-6008

*To place an order, please call: 800-851-9753, fax: 800-705-8443, or e-mail: [purchase@simpexmedical.com](mailto:purchase@simpexmedical.com)*