



## LARGE FRAGMENT

### Washers

Material: Stainless Steel

*Fabricated in Germany by Ortho Select® Implant Technology an ISO 13485:2003 Certified Company.*

Screw Sizes (mm)	Outer Diameter (mm)	Inner Diameter (mm)	Part No.
1.5 / 2.0	4.5	2.1	6099-0-0004
2.7 / 3.5	7.0	4.0	6099-0-0007
4.0 / 4.5	10.0	5.0	6099-0-0010
6.8	13.0	7.2	6099-0-0013
7.0	16.0	8.0	6099-0-0016