



LARGE FRAGMENT

Narrow Compression Plates 4.5 mm

Material: Stainless Steel
 Screw diameter: 4.5 mm
 Plate thickness: 4.0 mm
 Plate width: 12.0 mm
 Hole spacing: 16.0 mm - 25.0 mm
 Application: usually used for tibia shaft or small humerus shaft fixation.

Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.



Holes	Length (mm)	Part No.	Length (mm)	Minimal Contact Part No.
2	39	6012-0-6002	34	6012-0-4002
3	55	6012-0-6003	52	6012-0-4003
4	71	6012-0-6004	70	6012-0-4004
5	87	6012-0-6005	88	6012-0-4005
6	103	6012-0-6006	106	6012-0-4006
7	119	6012-0-6007	124	6012-0-4007
8	135	6012-0-6008	142	6012-0-4008
9	151	6012-0-6009	160	6012-0-4009
10	167	6012-0-6010	178	6012-0-4010
11	183	6012-0-6011	196	6012-0-4011
12	199	6012-0-6012	214	6012-0-4012
13	215	6012-0-6013		
14	231	6012-0-6014	250	6012-0-4014
15	247	6012-0-6015		
16	263	6012-0-6016	286	6012-0-4016
18	295	6012-0-6018		