



## LARGE FRAGMENT

### Broad Compression Plates 4.5 mm

Material: Stainless Steel  
Screw diameter: 4.5 mm  
Plate thickness: 4.8 mm  
Plate width: 16.0 mm  
Hole spacing: 16.0 mm - 25.0 mm  
Application: usually used for femur or humerus shaft fixation

*Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.*



Holes	Length (mm)	Part No.	Length (mm)	Minimal Contact Part No.
5	87	6013-0-6005		
6	103	6013-0-6006	106	6013-0-4006
7	119	6013-0-6007	124	6013-0-4007
8	135	6013-0-6008	142	6013-0-4008
9	151	6013-0-6009	160	6013-0-4009
10	167	6013-0-6010	178	6013-0-4010
11	183	6013-0-6011	196	6013-0-4011
12	199	6013-0-6012	214	6013-0-4012
13	215	6013-0-6013	232	6013-0-4013
14	231	6013-0-6014	250	6013-0-4014
15	247	6013-0-6015	268	6013-0-4015
16	263	6013-0-6016	286	6013-0-4016
17	279	6013-0-6017	304	6013-0-4017
18	295	6013-0-6018	322	6013-0-4018
20	327	6013-0-6020		