



LARGE FRAGMENT

Semi Tubular Plates

Material: Stainless Steel
 Screw diameter: 4.5 mm
 Plate thickness: 1.0 mm
 Plate width: 11.0 mm
 Hole spacing: 16.0 mm - xx.0 mm
 Application: usually used as additional plate on proximal or distal tibia.

Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.



Not to scale

Holes	Length (mm)	Part No.
2	39	6014-0-6002
3	55	6014-0-6003
4	71	6014-0-6004
5	87	6014-0-6005
6	103	6014-0-6006
7	119	6014-0-6007
8	135	6014-0-6008
9	151	6014-0-6009
10	167	6014-0-6010
11	183	6014-0-6011
12	199	6014-0-6012