



LARGE FRAGMENT

Condylar Buttress Plates

Material: Stainless Steel
Screw diameter: 4.5 mm
Plate thickness: 5.0 mm
Plate width: 16.0 mm
Hole spacing: 16.0 mm
Application: usually used for intra-articular distal femur fixation.

Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.



Holes	Length (mm)	Left Leg Part No.	Right Leg Part No.
7	158	6027-0-0007	6027-0-0008
9	190	6027-0-0009	6027-0-0010
11	222	6027-0-0011	6027-0-0012
13	254	6027-0-0013	6027-0-0014
15	286	6027-0-0015	6027-0-0016