



LARGE FRAGMENT

Large T-Plates

Material: Stainless Steel
Screw diameter: 4.5 mm
Plate thickness: 2.0 mm
Plate width: 17.0 mm
Application: usually used for proximal humerus or medial tibial plateau fixation.

Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.



Holes	Length (mm)	Part No.
3	68	6015-0-6003
4	84	6015-0-6004
5	100	6015-0-6005
6	116	6015-0-6006
7	134	6015-0-6007
8	148	6015-0-6008
9	164	6015-0-6009
10	180	6015-0-6010
12	212	6015-0-6012