



## LARGE FRAGMENT

### L - Buttress Plates

Material: Stainless Steel  
Screw diameter: 4.5 mm  
Plate thickness: 2.0 mm  
Plate width: 16.0 mm  
Application: usually used for lateral tibial plateau fixation.

*Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.*



Holes	Length (mm)	Left Leg Part No.	Right Leg Part No.
4	87	6017-0-6004	6018-0-6004
6	119	6017-0-6006	6018-0-6006
8	151	6017-0-6008	6018-0-6008
10	183	6017-0-6010	6018-0-6010