



LARGE FRAGMENT

Cobra Head Plates

Material: Stainless Steel
Screw diameter: 4.5 mm
Plate thickness: 5.0 mm
Hole spacing: 16.0 mm
Head width: 40.3 mm
Application: usually used for for hip fusion and arthodesis.

Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.



Holes	Length (mm)	Part No.
8	170	6001-0-6050
9	186	6001-0-6051
10	202	6001-0-6052
11	218	6001-0-6053