



LARGE FRAGMENT

Lateral Tibial Head Buttress Plates

Material: Stainless Steel
Screw diameter: 4.5 mm
Plate thickness: 4.0 mm
Plate width: 14.0 mm
Hole spacing: 16.0 mm
Application: usually used for lateral tibial plateau or proximal tibia fixation.

Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.



Holes	Length (mm)	Left Leg Part No.	Right Leg Part No.
5	118	6019-0-6005	6020-0-6005
7	149	6019-0-6007	6020-0-6007
9	181	6019-0-6009	6020-0-6009
11	214	6019-0-6011	6020-0-6011
13	246	6019-0-6013	6020-0-6013