



LARGE FRAGMENT

Proximal Tibial Plates

Material: Stainless Steel
Screw diameter: 4.5 mm
Plate thickness: 3.2 mm
Plate width: 16.0 mm
Hole spacing: 23.0 mm
Application: usually used for proximal tibial fixation.

Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.



Holes	Length (mm)	Left Leg Part No.	Right Leg Part No.
5	80	6029-0-0005	6029-0-0006
7	122	6029-0-0007	6029-0-0008
9	158	6029-0-0009	6029-0-0010
11	194	6029-0-0011	6029-0-0012
13	230	6029-0-0013	6029-0-0014