



LARGE FRAGMENT

Distal Tibial Plates

Material: Stainless Steel
Screw diameter: 4.5 mm
Plate thickness: 4.0 mm
Plate width: 14.0 mm
Hole spacing: 16.0 mm
Application: usually used for distal tibial fixation.

Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.



Holes	Length (mm)	Left Leg Part No.	Right Leg Part No.
5	80	6030-0-0005	6030-0-0006
7	122	6030-0-0007	6030-0-0008
9	158	6030-0-0009	6030-0-0010
11	194	6030-0-0011	6030-0-0012
13	230	6030-0-0013	6030-0-0014