



LARGE FRAGMENT

Distal Femoral Plates

Material: Stainless Steel
Screw diameter: 4.5 mm
Plate thickness: 3.2 mm
Plate width: 16.0 mm
Hole spacing: 24.0 mm
Application: usually used for distal femoral fixation.

Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.



Holes	Length (mm)	Left Leg Part No.	Right Leg Part No.
5	80	6028-0-0005	6028-0-0006
7	122	6028-0-0007	6028-0-0008
9	158	6028-0-0009	6028-0-0010
11	194	6028-0-0011	6028-0-0012
13	230	6028-0-0013	6028-0-0014