



## LARGE FRAGMENT

### Lengthening Plates 8 Holes

Material: Stainless Steel  
 Screw diameter: 4.5 mm  
 Application: usually used for fixation after prolongation osteotomy.

*Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.*



Lengthening by (mm)	Length (mm)	Narrow Part No.	Broad Part No.
Plate thickness:		3.6 mm	4.5mm
Plate width:		12.0 mm	16.0 mm
Hole spacing:		12.0 mm	12.0 mm
30	135	6033-0-7135	6012-0-7135
40	145	6033-0-7145	6012-0-7145
50	155	6033-0-7155	6012-0-7155
60	165	6033-0-7165	6012-0-7165
70	175	6033-0-7175	6012-0-7175
80	185	6033-0-7185	