



## LARGE FRAGMENT

### Lengthening Plates 10 Holes

Material: Stainless Steel  
Screw diameter: 4.5 mm  
Plate thickness: 4.5 mm  
Plate width: 16.0 mm  
Hole spacing: 12.0 mm  
Application: usually used for fixation after prolongation osteotomy.

*Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.*



Lengthening by (mm)	Length (mm)	Broad Part No.
50	180	6013-0-7180
60	190	6013-0-7190
70	200	6013-0-7200
80	210	6013-0-7210
90	220	6013-0-7220
100	230	6013-0-7230
110	240	6013-0-7240
120	250	6013-0-7250