

LARGE FRAGMENT


Cortex Screws 4.5 mm

Material:	Stainless Steel
Thread Diameter:	4.5 mm
Core Diameter:	3.0 mm
Head Diameter:	8.0 mm
Hex Diameter:	3.5 mm
Thread Pitch:	1.75 mm
Pilot Drill:	3.2 mm
Gliding Hole Drill:	4.5 mm
Tap:	4.5 mm
Screwdriver Hex:	3.5 mm
Application:	usually used for interfragmentary compression according to the lag-technique, for plate fixation, or for both.

Fabricated in Germany by Ortho Select[®] GmbH Implant Technology; ISO 13485:2003 Certified.




Length (mm)	4.5 mm Standard Part No.	4.5 mm Self-Tapping Part No.
14	6080-0-6014	6081-0-6014
16	6080-0-6016	6081-0-6016
18	6080-0-6018	6081-0-6018
20	6080-0-6020	6081-0-6020
22	6080-0-6022	6081-0-6022
24	6080-0-6024	6081-0-6024
26	6080-0-6026	6081-0-6026
28	6080-0-6028	6081-0-6028
30	6080-0-6030	6081-0-6030
32	6080-0-6032	6081-0-6032
34	6080-0-6034	6081-0-6034
36	6080-0-6036	6081-0-6036
38	6080-0-6038	6081-0-6038
40	6080-0-6040	6081-0-6040



LARGE FRAGMENT

Cortex Screws 4.5 mm



Length (mm)	4.5 mm Standard Part No.	4.5 mm Self-Tapping Part No.
42	6080-0-6042	6081-0-6042
44	6080-0-6044	6081-0-6044
46	6080-0-6046	6081-0-6046
48	6080-0-6048	6081-0-4048
50	6080-0-6050	6081-0-6050
52	6080-0-6052	6081-0-6052
54	6080-0-6054	6081-0-6054
56	6080-0-6056	6081-0-6056
58	6080-0-6058	6081-0-6058
60	6080-0-6060	6081-0-6060
65	6080-0-6065	6081-0-6065
70	6080-0-6070	6081-0-6070
75	6080-0-6075	6081-0-6075
80	6080-0-6080	6081-0-6080
85	6080-0-6085	6081-0-6085
90	6080-0-6090	6081-0-6090
95	6080-0-6095	6081-0-6095
100	6080-0-6100	6081-0-6100
105	6080-0-6105	6081-0-6105
110	6080-0-6110	6081-0-6110
115	6080-0-6115	6081-0-6115
120	6080-0-6120	6081-0-6120
125	6080-0-6125	6081-0-6125
130	6080-0-6130	6081-0-6130
135	6080-0-6135	6081-0-6135
140	6080-0-6140	6081-0-6140