

LARGE FRAGMENT

Cancellous Screws 6.5 mm

Material:	Stainless Steel
Thread Diameter:	6.5 mm
Core Diameter:	3.0 mm
Shaft Diameter:	4.5 mm
Head Diameter:	8.0 mm
Hex Diameter:	3.5 mm
Thread Pitch:	2.75 mm
Pilot Drill:	3.2 mm
Gliding Hole Drill:	4.5 mm
Tap:	6.5 mm
Screwdriver Hex:	3.5 mm
Application:	The partially threaded screws are usually used as lag screws in the epiphyseal and metaphyseal areas, where the cortex is thin and the bone cancellous. The two thread lengths (16 mm, 32 mm) are available in order to accommodate different sizes of the bone fragment.

The fully threaded screws are usually used with plates.

Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.



Length (mm)	16 mm Thread Length Part No.	32 mm Thread Length Part No.	Fully Threaded Part No.
20			6082-0-8020
25	6082-0-6025		6082-0-8025
30	6082-0-6030		6082-0-8030
35	6082-0-6035		6082-0-8035
40	6082-0-6040	6082-0-7040	6082-0-8040
45	6082-0-6045	6082-0-7045	6082-0-8045
50	6082-0-6050	6082-0-7050	6082-0-8050
55	6082-0-6055	6082-0-7055	6082-0-8055

LARGE FRAGMENT

Cancellous Screws 6.5 mm



Length (mm)	16 mm Thread Length Part No.	32 mm Thread Length Part No.	Fully Threaded Part No.
60	6082-0-6060	6082-0-7060	6082-0-8060
65	6082-0-6065	6082-0-7065	6082-0-8065
70	6082-0-6070	6082-0-7070	6082-0-8070
75	6082-0-6075	6082-0-7075	6082-0-8075
80	6082-0-6080	6082-0-7080	6082-0-8080
85	6082-0-6085	6082-0-7085	6082-0-8085
90	6082-0-6090	6082-0-7090	6082-0-8090
95	6082-0-6095	6082-0-7095	6082-0-8095
100	6082-0-6100	6082-0-7100	6082-0-8100
105	6082-0-6105	6082-0-7105	6082-0-8105
110	6082-0-6110	6082-0-7110	6082-0-8110
115	6082-0-6115	6082-0-7115	6082-0-8115
120	6082-0-6120	6082-0-7120	6082-0-8120
125	6082-0-6125	6082-0-7125	6082-0-8125
130	6082-0-6130	6082-0-7130	6082-0-8130
135	6082-0-6135	6082-0-7135	
140	6082-0-6140	6082-0-7140	
145		6082-0-7145	
150		6082-0-7150	