



LARGE FRAGMENT

Malleolar Screws 4.5 mm

Material:	Stainless Steel
Thread diameter:	4.5 mm
Core diameter:	3.0 mm
Head diameter:	8.0 mm
Hex diameter:	3.5 mm
Thread pitch:	1.25 mm
Pilot drill:	3.2 mm
Screwdriver hex:	3.5 mm
Application:	usually used to fix fragments in the metaphysis area, malleolus, distal humerus, and lesser trochanter.

Fabricated in Germany by Ortho Select® GmbH Implant Technology; ISO 13485:2003 Certified.



Length (mm)	Part No.
20	6080-0-5020
25	6080-0-5025
30	6080-0-5030
35	6080-0-5035
40	6080-0-5040
45	6080-0-5045
50	6080-0-5050
55	6080-0-5055
60	6080-0-5060
65	6080-0-5065
70	6080-0-5070
75	6080-0-5075
80	6080-0-5080